



Easily tell if the valve is open or closed with our patented visual indicator feature.

## OXY-TRU<sup>®</sup> MEDICAL GAS VALVE


Give your cylinder fill plant personnel a quick visual confirmation on the valve's position (closed or open) by simply looking down the fill manifold line.

### SPECIFICATIONS

- Inlet Straight Thread – ¾ - 16 UNF straight thread
- Outlet – CGA 870 USP Oxygen
- Cleaned for Medical Oxygen service
- Valve Service Pressure = 2,015 psi
- Chrome plated brass
- Design burst pressure = 6,750 psi
- Meets CGA V-9

### SAFETY INFO

- Rupture Disc – 11083E-12C 3,360 psi
- Flow Coefficient (CV) = 0.14



### One Resource. Zero Headaches.<sup>®</sup>

Your one resource for gas cylinders, cryogenic vessels, small & large bulk tanks, valves, and cryo repair services.

*(Customers outside U.S. call +1-630-844-8800.)*

U.S. TOLL FREE **888-429-5832**

[www.cyl-tec.com](http://www.cyl-tec.com)  **Cyl-Tec**<sup>®</sup>